

Glenelg J-48
Sample Overview

SAMPLES

Sample Type	Interval (m)	# of Samples
Washed Cuttings	380 – 5,250	900
Unwashed Cuttings	380 – 5,250	906
Sidewall Core	468 – 5,107	464
Canned Cuttings (dried)	380 – 5,240	507

Slides:	Interval (m)	# of Slides	Sample Source
Micropaleo	370 – 5,250	167	cuttings
Micropaleo	468 – 4,750	218	sidewall core
Palynology	370 – 5,250	167	cuttings
Palynology	5,086 – 5,092	4	company cuttings

Core:	Interval (m)	Recovery
Core #1	5,179.0 – 5,188.0	9.9
Core #2 (J-48A)	5,085.0 – 5,095.0	8.86

Fluids:	Interval (m)	Recovery	Recovered from
Test #			
DST #8, Zone 7	3,491 – 3,495.5m	condensate	separator
DST# 9, Zone 8	3,061.9 - 3,064.9m	condensate	separator